



Assessing Aerosol Exposure from High Speed Cell Sorters

Jeff Potts, MPH

Occupational Health & Safety Specialist

National Institutes of Health

Division of Occupational Health and Safety

Department of Health and Human
Services

Introduction

- Cytometry is the measurement of physical and/or chemical characteristics of cells.
- Flow Cytometry is a process in which such measurements are made while the cells pass, single file, through the measuring apparatus in a fluid stream.
- Measure the presence or absence of cell surface receptors, intracellular proteins or nucleic acids, and functional changes
 - Performed rapidly (up to 30,000 cells/sec)
 - Performed quantitatively
- Cell Sorting
 - Isolate cell subsets for subsequent functional or genomic/proteomic analysis

Introduction

- Droplet Sorting
 - Advantages:
 - Speed (30,000 cells/sec)
 - Up to 4 populations simultaneously sorted
 - Sorted cells relatively concentrated
 - Disadvantages
 - Aerosols generated
 - Expensive
 - High speed (high pressure) sorting may damage cells

How does a droplet cell sorter sort?

- Physical fluidics principle:
 - A stream of fluid is hydrodynamically unstable
 - It will break up into droplets that have a smaller surface area and lower surface tension
 - Droplet formation pattern occurs in an unpredictable manner


Droplet formation in sink



Introduction

- Droplet generation
 - Droplet pattern can be stabilized
 - Nozzle vibration induces an acoustic wave in stream
 - Droplets form at troughs of waves
 - Breakoff therefore occurs in the same place at a given frequency of vibration

Droplet Controlled by Nozzle Vibration Frequency





Introduction

- This preliminary study was conducted to ascertain the feasibility of using particle sampling to detect aerosolization created by cell sorting activities.
 - If particle sampling is feasible it could be used to:
 - Quickly assess AMS prior to infectious cell sorting.
 - Identify potential hazardous operating procedures
 - Evaluate cytometer design and modifications

Introduction

- Currently two methods for measuring aerosol generation are highlighted in the International Society for Analytical Cytology (ISAC) Biosafety Standard for Sorting of Unfixed Cells published in January 2007.
 - T4 bacteriophage method
 - Fluorescent bead method

Introduction

- T4 bacteriophage method
 - T4 bacteriophage solution is ran through the flow cytometer.
 - T4 bacteriophage susceptible *Escherichia coli* lawn plates are distributed through out the work area.
 - Plates are incubated for 24 hrs then read for plaques.
 - Number of plaques correlate with the T4 bacteriophage aerosol deposition.

Introduction

- Fluorescent bead method
 - Fluorescent beads are suspended in a phosphate buffered saline solution, then ran through the flow cytometer.
 - Microscope slides are places through out the work area.
 - Slides are then examined for aerosolization under a fluorescent microscope.
 - Number of identified fluorescent particles identified correlated to aerosol generation.

Introduction

- Concerns Associated with current methods.
 - T4 bacteriophage method
 - 24 incubation required and extensive prep time required
 - Knowledge of microbiology required
 - Contamination of flow cytometer
 - Fluorescent bead method
 - Individual variability
 - Time consuming
 - Limited particle size
 - Diminished fluorescence

Introduction

- Aerodynamic Particle Sampling
 - Particle sampling not new
 - Historical concerns associated with particle sampling
 - Interference from background particles
 - Inability to differentiate particle size
 - Particle impaction in sampling tube
 - Conservative due to assumption that particles represent biological material

Introduction


- Aerodynamic Particle Sampling
 - Particles are pulled into instrument and characterized by light scattering technology.
 - Particles are counted and categorized by aerodynamic size.
 - Instrument categorizes sizes between 0.5 microns (μm) and 20 μm .
 - Samples collected and averaged during 20 second intervals.

Methodology




Methodology

- Location 0 – Stream source
- Location 1 – 7 moving away from source in 1.3 cm increments
- Location 8-11 increments increased to 6.1 cm




Results



Graph 1: Room No operation


Results



Graph 2: Inside Biological Safety Cabinet No Operation




Results




Graph 3: Particle Counts Moving Away from Stream During Replicated Failure

Results




Graph 4: Particles @ Location 0 During Replicated Failure

Results




Graph 5: Particles @ Location 11 During Replicated Failure

Results




Graph 6: Particles @ Location 0 During Replicated Failure @ 70 psi

Results




Graph 7: Particles @ Location 0 During Replicated Failure @ 20 psi

Results



Graph 8: Measurements w/o Stream Disruption

Results



Graph 9: Measurements w/ Door Closed

Discussion

- When background particles are eliminated aerodynamic particle sampling is capable of assessing potential inhalation exposure from cell sorting.
 - Further studies will be conducted utilizing a particle sampler which has capabilities of detecting particles tagged with a fluorescent compound.
 - An ultraviolet (UV) laser will be used to identify particles that contain a UV-excitable dye (Clear Blue), which has been dissolved in the fluidics system of the cell sorter.
 - This will allow counting and sizing of aerosols unambiguously identified as being generated by the cell sorter versus those in the ambient air.

Discussion

- Goals of Future Studies
 - Evaluate current methods for test AMS
 - Develop safe operating procedures in the event of stream disruption
 - Primary Barriers
 - Assess equipment modifications for potential aerosolization hazards
 - Provide accurate quantification of potential aerosolization hazards to be used by safety specialist in conducted risk assessments.
 - Personal Protective Equipment

Discussion

- Problems associated with using microbial exposure to determine risk of developing disease.
 - Particles generated may or may not carry infectious agent
 - Fixed samples as an alternative
 - Individual susceptibility
 - Immunizations
 - Virulence of infectious particles



Acknowledgments

LT Gino K. Begluitti, MPH, CHMM
Industrial Hygienist
National Institutes of Health
Division of Occupational Health and Safety

Kevin Holmes, PhD
Chief, Flow Cytometry Operations
National Institutes of Health
National Institute of Allergies and Infectious Diseases

Angela Ragan
Fellow
National Biosafety & Biocontainment Training Program
National Institutes of Health